




### Using a Platinum stairlift is very simple.




Sit on the revolutionary Ergo chair that can be fully adjusted to ensure the perfect fit for you at installation.



At the top you can swivel the chair onto the landing to allow safe transfer on and off the stairlift.




To operate the stairlift simply use the ergonomic joystick in the end of the arm by lightly pressing the joystick in the direction you wish to travel.



When not in use the stairlift can be parked away on the stairs to leave space for other family members.



# Discover the unique features.



#### **Ergonomic joystick**


Can be used as a button or lever to make the stairlift easy to operate in a natural and comfortable position even with limited mobility.





#### Security key switch

under the armrest, to prevent unauthorized use of your stairlift.



# Independent lever linked footrest raiser

to open and close the footrest easily when you need it without having to bend down (powered option also available). Also allows easy transfer for wheelchair users.



#### **ERGO** plus seat

designed to optimize the position of the legs when seated.




# Digital diagnostic display

built into the armrest, to check all the information and diagnostics about you stairlift in real time.



### Retractable reel seat belt

to provide additional security to move safely.



#### Seat swivel lever

To turn you at the top of the stairs simply manually lift the lever up. Powered swivel option also available which works off the ergonomic joystick.



These would immediately stop the stairlift in the event of any obstructions on the stairs during travel.









### and easy to use.

- Start delay to prevent accidental operation
- Choice of colours. Change the upholstery and look new in seconds.
- Battery operated so works in the event of a power cut
- User weight capacity up to 160 kg on curve and 140 kg on straight.

# Upholstery colors. Quality fabrics and tailored finishes.










**Grey** Vinyl

Red Vinyl

Blue Vinyl

Beige Vinyl



## Curved or straight,

### 




Power supply	Batteries, 2 x 12 volt. consumption 0,014 kWh
Recharge	Automatic, in the charging points
Speed	0,12 m/sec.
Wiring	Internal with charging plates at the ends of the rail
Manufacture	Aluminium single rail
Anchorage	Mounting step through telescopic brackets. Wall mounting optional.
Traction	Rack and pinion with OSG control
Diagnostic	built into armrest, internal memory and operating report
Rail slope	25 - 62° max
Engine	Battery supply, 0,26kw
Recharge	External unit, 24 volt in two stages
Warranty	24 months
Capacity	140 kg

(footplate adjustment)



# the Chair is always ERGO.



The measurements, in millimeters, are nominal values.
May vary depending on the configuration of the stair.
The instructions for use are supplied with each system. Platinum and Ergo trademarks and logos are registered.

